

7 H F K Q L F D O ' D W D S P K D W H W R U 3 O D W H P H

### Application:

The capmats for 3ate medio enables to cover the round well plates tight and impermeable up to the temperature of -25°C. It has been developed in line with the regulations of the Society for Biomolecular Screening (SBS).

### Features:

- x Capmats for 3ate medio is produced out of ethylene vinyl acetate.
- x Each capmat is equipped with an alphanumeric code for easy and fast identification of each sample.

### Article no.:

### Dimensions:

length: 120 mm (+/-0,5)  
width: 78 mm (+/-0,5) 4,5  
height: mm (+/-0,5)


**Wells:** 6 x 6 mm, centered with a distance of 9 mm (SBS-format)

**Edges:** round edges with an extension on the right side



Cod: 243114  
capmat for Plate medio

**Dettaglio**  
Cod: 243114  
capmat for Plate medio



Specifications are subject to change without prior notice.  
In case of any questions, please contact our sales representative via phone, e-mail or our website.  
[www.picov.eu](http://www.picov.eu)

Manufactured in Germany by Ritter GmbH - Disposable capmat for Plate medio